



Model 330 Series
Users Manual

WEEE/RoHS Compliance Statement


EU Directives WEEE and RoHS

To Our Valued Customers:

We are committed to being a good corporate citizen. As part of that commitment, we strive to maintain an environmentally conscious manufacturing operation. The European Union (EU) has enacted two Directives, the first on product recycling (Waste Electrical and Electronic Equipment, WEEE) and the second limiting the use of certain substances (Restriction on the use of Hazardous Substances, RoHS). Over time, these Directives will be implemented in the national laws of each EU Member State.

Once the final national regulations have been put into place, recycling will be offered for our products which are within the scope of the WEEE Directive. Products falling under the scope of the WEEE Directive available for sale after August 13, 2005 will be identified with a "wheelie bin" symbol.

Two Categories of products covered by the WEEE Directive are currently exempt from the RoHS Directive - Category 8, medical devices (with the exception of implanted or infected products) and Category 9, monitoring and control instruments. Most of our products fall into either Category 8 or 9 and are currently exempt from the RoHS Directive. We will continue to monitor the application of the RoHS Directive to its products and will comply with any changes as they apply.



- **Do Not Dispose Product with Municipal Waste**
 - **Special Collection/Disposal Required**

KDS Model 330 Series

Contents

Specifications.....	1
Safety	2
Features.....	3-4
Operating Instructions	
Power Switch	5
Syringe Loading.....	5
Menu Features	5
Keypad Functions.....	6
Fast Forward, Fast Reverse	6
Syringe Diameter Entry	6
Mode Selection	7
Volume Setting	7
Rate Setting.....	8
Power Up (run or stop)	8
Run/Stop.....	8
Change or Review Volume Setting While Running	8
Clearing a Stall Condition	9
NV Ram Failure	9
Fuses	9
Maintenance	10
Standard Syringe Diameters	11
Japanese Syringe Diameters	12
Standard Minimum and Maximum Flowrates.....	13
Japanese Minimum and Maximum Flowrates.....	14
Limited Warranty.....	15

KDS Model 330 Series

Specifications

Model	Model 330
Syringe size	10 ml Popper 5127
Electrical Rating	US model 115 V~, 0.25 A CE model 230 V~, 0.16 A
Fuse	5 x 20 mm, 250 V~, Slow Blow, 0.25 A
Voltage operating range	US model 110-120 V~, 50/60 Hz CE model 200-240 V~, 50/60 Hz
Drive mechanism	Microprocessor controlled stepper motor 1/2 - 1/16 microstepping, driving a leadscrew through a belt and pulley drive mechanism
Pusher advance per microstep	(1/16 step) - 0.165 micron (or .0000064 inch)
Volume per microstep	(1/16) with 60 ml BD syringe - 0.0919 μ l
Max stepping rate	1600 (1/2 step) /sec
Min step rate	1 (1/16 step) /120 secs
Linear Travel Rate	Min: 4.95×10^{-4} cm/hr Max: 12.67 cm/min
Flowrate range	2.757 μ l/hr to 70.56 ml/min (60 ml syringe)
Nominal linear force	35 lb
Dimensions	23 x 7.5 x 5.5 in (58 x 19 x 14 cm)
Weight	16 lb (7.5 kg)
Atmospheric Specifications	
Temperature	5°C - 40°C (41°F - 104°F)
Humidity	20% - 80% RH
Mode of Operation	Continuous
Classification	Class I
Pollution Degree	2
Installation Category	II
Output	N/A
Physiological Effects	N/A
Cooling Conditions	No special considerations
Mechanical Stability	No special considerations
Protective Packaging	No special considerations
Earth Terminals	No External connections required
Removable Protective Means	N/A
Supplier Name	KD Scientific Inc.
Address	84 October Hill Road, Holliston, MA 01746-1388

KDS Model 330 Series

General Safety Summary

Please read the following safety precautions to ensure proper use of your syringe pump. To avoid potential hazards and product damage, use this product only as instructed in this manual. If the equipment is used in a manner not specified by the manufacturer, the protection provided by the equipment may be impaired.

To Prevent Hazard or Injury:

USE PROPER POWER SUPPLY

The pump is supplied with an approved power supply and line cord.

USE PROPER LINE CORD

Use only the line cord shipped with the product and make sure line cord is certified for country of use.

GROUND PRODUCT

This product is grounded through the grounding conductor of the power cord. To avoid electric shock, use only approved line cord with the product and insure it is connected to earth ground.

MAKE PROPER CONNECTIONS

Make sure all connections are made properly and securely.

ORIENT EQUIPMENT PROPERLY

Do not position the equipment such that it is difficult to reach the disconnecting device.

OBSERVE ALL TERMINAL RATINGS

Review the operating manual to learn the ratings on all connections.

AVOID EXPOSED CIRCUITRY

Do not touch any electric circuitry inside the product.

DO NOT OPERATE WITH SUSPECTED FAILURES

If damage is suspected on or to the product do not operate the product. Contact qualified service personnel to perform inspection.

AVOID PINCH HAZARD

A pinch hazard may exist between the Pusher Block and End Blocks. Avoid placing fingers between these points while the pump is running.

OBSERVE ALL WARNING LABELS ON PRODUCT

Read all labels on product to ensure proper usage.



CAUTION
Refer to Manual



*Protective Ground
Terminal*

CAUTION

This pump is not registered with the FDA and is not for clinical use on human or veterinary patients. It is intended for research use only.

CAUTION
FOR RESEARCH USE ONLY
NOT FOR CLINICAL
USE ON PATIENTS

KDS Model 330 Series

Features

The model KDS330 is designed for emulsifying proteins in viscous fluids. This pump retains the software features of the infusion/withdrawal syringe pump, KDS210. For emulsion making, the Continuous mode is selected and the pump is specifically designed for a 10 ml glass syringe and **emulsion needle**.

Operation of the pumps is simplified by using a keypad to select features from a menu on the alphanumeric LCD display.

All control functions are performed automatically by the pump microcontroller and are based on the syringe diameter and linear motion of the pusher block to deliver a known volume. After entering the syringe diameter, either directly or from a table in memory, a dispense volume and flow rate can be entered, and then all calibration and control functions are performed by the pump automatically.

1) Syringe identification

Look up Table

The pump contains a table of standard syringes arranged by manufacturer, material and size. Once the syringe is identified in the table the pump automatically enters the appropriate diameter.

Direct Entry

If the syringe used is not included in the table, the internal diameter of the syringe barrel can be measured in millimeters and entered directly from the keypad.

2) Mode of operation - Continuous operation

This mode is used in emulsion making. The pump cycles from infusion to withdrawal continuously. The volume is identical in infuse and withdrawal directions (the modes listed below are available but are not needed for emulsion making).

Infusion

Rate and volume settings: pump infuses to the set volume and stops.

Rate setting only: pump runs until manually stopped or stalls.

Withdrawal

Rate and volume settings similar to above.

Infusion/Withdrawal

Infusion automatically followed by withdrawal. Rate and volume settings can be made independently for infusion and withdrawal.

Withdrawal/Infusion

Withdrawal immediately followed by infusion. Separate settings for rate and volume can be made for withdrawal and infusion.

KDS Model 330 Series

Features (Continued)

3) Volume

In "Continuous" mode a target volume must be entered. The pump displays an initial volume of zero and increases as the dispense proceeds to the target volume. The target volume can be reviewed or changed as the pump continues to operate

[In other modes the infusion and withdrawal volumes can be set independently.]

4) Flow Rate

In "Continuous" mode only one flowrate is entered and the pump cycles back and forth at this rate.

[In other available modes, the infusion rate and, where applicable, the withdrawal rate can be set independently and can be changed while the pump is running. After the operating mode selection is made the program will prompt only for the relevant rates associated with that mode.]

5) Stall Detection

The motor is monitored by an optical encoder to confirm the programmed movement. If the back pressure increases due to jamming or flow restriction then the motor may stall. Stall detection by the encoder results in a pump shutdown. The display will read "Stalled". The Stall message can be cleared with the Select key.

6) Power Disruption

When power is returned after a temporary power disruption the pump can be programmed to resume operation or remain stopped. However, if a dispense volume is set then the pump always remains stopped.

7) Non-volatile Memory

All operational settings are stored in non-volatile memory for convenience, and are used to set the pump when first switched on.

8) Selection of Rate and Volume Units

Units of volume (μ l or ml) and flowrate (μ l/ml per min/hr) can be changed if required.


KDS Model 330 Series

Operating Instructions

POWER SWITCH

The power switch is located on the left corner of the rear panel.

SYRINGE LOADING



To simplify syringe loading the pusher block can be disengaged from the leadscrew by pressing the brass button and manually moved along the guide rods. Alternately, the Fast forward, Fast reverse feature can be used (press Run and Arrow key simultaneously).

One syringe is filled only slightly to remove any air in the "dead stop" volume of the syringe. The second syringe is filled with fluid mixture and any air is expelled. The emulsion needle is attached to this second syringe and any air bubbles are expelled from the emulsion needle. The first syringe is then attached to the opposite opening of the emulsion needle, creating a closed system with the two syringes apposed, and containing fluid devoid of air pockets. The two clamps on the syringe blocks should be rotated and the syringe barrels placed in the "V" of the syringe holders. With the syringes in place the retaining clamps can be rotated and tightened down on the barrel to hold it securely in place.

Release the pusher block mechanism and slide the block along the guide rods until the syringe plunger fits securely in the pusher block.

MENU FEATURES

When the pump is first turned on the LCD will display the most recent Rate or Volume setting. Pressing **select** repeatedly will always return the display to the main menu.

The main menu consists of seven variables, three of which are displayed at all times. The center variable pulses to indicate that this option can be reviewed or changed.

The menu acts as a continuous loop and the **arrow keys**, **←** **→** are used to move around the loop. Press the **select** key to select the variable to be changed.

← → DIA ← → TABLE ← → RATE ← → VOL ← → MODE ← → RS232C ← → PWR UP ← →

KDS Model 330 Series

Operating Instructions (Continued)

KEYPAD FUNCTIONS

A keypad is used for selection of features from the displayed menu and numerical entries.

KEYPAD



These keys move the displayed menu left and right respectively.

select

Used to select the highlighted menu feature (flashing).

enter

The numerical parameters of the feature selected are set, then enter places all the settings in memory ready for operation. All settings are stored in non-volatile memory and will be used next time the pump is operated.

0 - 9, .

Numerical and decimal entry keys.

run/stop

Starts and stops the motor. Acts as a pause during a dispense.

FAST FORWARD, FAST REVERSE

Keys run and \rightarrow or \leftarrow pressed simultaneously causes the pump to run at maximum speed. This feature can be used for loading, purging and reversing out of a stall condition.

SYRINGE DIAMETER ENTRY

The pump is designed for use with Popper & Sons 10ml glass syringes with metal Luer lock tip (p/n 5128) in conjunction with a micro-emulsifying needle, 18g x 2.875" (p/n 7973).

The pump must be calibrated by identifying the internal diameter of the syringe used. Once entered this data is stored in non-volatile memory and need be modified only when a different syringe is used.

Table

Use the **arrow** keys to move TABLE to the center of the display and press select to **select** this feature. The display now reads abbreviated names of syringe manufacturers and the type of syringe (plastic or glass).

Use the **arrow** and **select** keys to scroll through the list and select the manufacturers name of the syringe used.

The display now lists the size of syringe. Again, use the **arrow** and **select** keys to identify and enter the syringe used.

Diameter

If the syringe used is not listed in the table of syringes then the internal diameter of the syringe must be measured and entered directly. Scroll through the main menu and select **DIA**.

The display now prompts for entry of the syringe barrel internal diameter measured in **millimeters**. Enter by using the numerical keys and **enter**.

Note: If the diameter is changed the volume and rate settings are set to zero.

KDS Model 330 Series

Operating Instructions (Continued)

MODE SELECTION

Select MODE from the main menu and then scroll through the options displayed and select the mode required – **Continuous**

Continuous

The pump first infuses and then withdraws, and then cycles continuously. Only one volume setting for infusion and withdrawal is permitted.

The menu now prompts for volume and rate settings relevant to the mode selected.

VOLUME SETTING

Volume can be reached directly from the main menu or will be prompted after mode selection.

Volume prompt: Vol.: 00.00ml >

- 1) Enter target volume from the numerical keypad.
- 2) If the units and value displayed are correct, **enter**.

The underlined display, or pulsing display segment, indicates that this parameter can be changed. The arrow symbol signifies direction of travel (< indicates withdrawal) and pulses when the pump is running.

Note: When dispensing, the displayed volume increments in units of the last significant figure of the volume entered. Therefore, to increment by .01 the volume should be set at 1.00. If the volume is set at 1.000 then the volume increments in units .001.

Note: Volume setting = 0.0

- a) Infusion or Withdrawal modes.
Volume = 0 is interpreted as no volume and the pump will run until manually stopped or a stall occurs. The LCD will display the rate setting.
- b) In continuous mode a volume setting must be entered. If no volume is entered the display returns to the mode selection menu to select the appropriate mode.
- c) Unit setting

prompt: Vol.: 00.00 ml <

The units displayed can be changed if required.

- 1) Use the → arrow key to move the pulsing indicator to the units displayed.
- 2) Continue to use the → arrow key to scroll through the possible units.
- 3) The ← arrow key will move the active display back to the numerical value.
- 4) When the correct value and units are displayed press **enter**.

Possible units are: µl and ml microliter and milliliter

KDS Model 330 Series

Operating Instructions (Continued)

RATE SETTING

Display reads: Rate 00.01 m³/h >

- 1) Enter the flowrate value required with the numerical keypad
- 2) If necessary, change the units using the → key to move to and scroll through the possible units.
- 3) When the displayed settings are correct press **enter**.

Note: If the number entered exceeds the maximum flow rate possible then the pump displays the maximum feasible rate. To continue enter a rate smaller than the maximum.

Note: To check the maximum possible rate enter 9's to the required decimal position. For example, enter 99.9 and the maximum displayed is 12.3 whereas, if 99.99 is entered, then the maximum displayed is 12.34.

POWER UP run or stop

This option is only applicable when **no dispense volume** is selected. When power returns after an interruption the pump can resume operation (select RUN) or remain stopped (select STOP).

If the pump resumes operation the rate display will flash to indicate that a power interruption has occurred. Press **select** to clear the display to resume normal operation.

RUN/STOP

After all settings are made the pump can be started or stopped by a single press of the **run/stop** key. During a volume dispense the **stop** acts as a "pause" and **run** will resume the dispense.

CHANGE OR REVIEW VOLUME SETTING WHILE RUNNING

While the pump continues to run press **select** to return to the main menu. Scroll through the menu and select Volume to display the set dispense volume.

a) No volume change

Press **select**. The display returns to the incrementing display volume.

b) Volume change

- 1) Make the changes with the numerical keypad and **enter**.
- 2) The display moves to RATE, permitting a change if required. Use the numerical keypad and enter to make changes (or leave unchanged). The pump immediately changes to the new flowrate, if changed, and the volume continues to increment, uninterrupted by the review process, to the target dispense volume, and will stop automatically when completed.

Note: If the volume is changed to a volume smaller than the volume already accumulated then the pump will stop as soon as the new, smaller target volume is entered.

KDS Model 330 Series

Operating Instructions (Continued)

CLEARING A STALL CONDITION

Should a stall occur the pump motor is stopped to prevent damage.

To clear the display press **select**.

To move the stalled mechanism use the **fast forward** or **fast reverse** to move the pusher block. Using the fast forward or fast reverse feature is not only the most simple way to deal with the stall, it also reduces potential damage to the cam mechanism which releases the halfnut from the leadscrew. Alternatively, depress button on pusher block.

NV RAM FAILURE

If the settings in the non-volatile memory become corrupted the display will read "NV Ram Failure" and the pump will not operate.

To recover from this condition the pump must be powered down and then turned on again. The pump will be re-initialized to the default settings and can now be programmed as normal.

The "NV Ram Failure" message can also be cleared by pressing select and programming a new flow rate. The pump should then be turned off and on to save the settings.


Fuses

The fuses are located, in the power entry module on the rear panel. The linecord must be removed first to gain access to the fuse holder.

Fuse 5 x 20 mm, 250V~ slow blow, 0.25 A

Voltage Selector (CE version only)

If it is necessary to change the input voltage selection, disconnect the line cord from the entry module on the rear panel. Use a flat bladed screwdriver to open the Fuse Holder access door. Remove the Fuse Holder, flip over, and reinstall. Close the access door. The new input voltage selection should be visible through the door window. Intall a proper line cord certified for the country of use.



KDS Model 330 Series

Maintenance

Maintenance is required only for the moving mechanical parts, which should be kept clean and lubricated. Occasionally, a small amount of light machine oil should be applied to the guide rods and a small amount of grease or oil to the leadscrew.

Solvents of any type should never be used to clean the pump. A mild detergent solution may be used to clean the keypad.

KDS Model 330 Series

Standard Table of Syringe Diameters

(1) "Air-Tite "All Plastic		(6) Ranfac	
1cc	4.70mm	2cc	9.12mm
2.5	9.70	5	12.34
5.0	12.48	10	14.55
10	15.89	20	19.86
20	20.00	30	23.20
30	22.50	50	27.60
50	28.90		
(2) Becton Dickinson		(7) Scientific Glass Engineering	
Interim, WW design, Plastipak		SGE	
1cc	4.70mm	25 μ l	0.73mm
3	8.59	50	1.03
5	11.99	100	1.46
10	14.48	250	2.30
20	19.05	500	3.26
30	21.59	1ml	4.61mm
60	26.60	2.5	7.28
		5	10.30
		10	14.57
(3) Becton Dickson		(8) Sherwood - Monojet Plastic	
Glass - all types		1cc	4.65mm
0.5cc	4.64mm	3	8.94
1	4.64	6	12.70
2.5	8.66	12	15.90
5	11.86	20	20.40
10	14.34	35	23.80
20	19.13	50	26.60
30	22.70		
60	28.60	(9) Terumo	
		1cc	4.73mm
(4) Hamilton		3	9.00
1000-Series Gastight		5	13.04
10 μ l	0.46mm	10	15.79
25	0.73	20	20.18
50	1.03	30	23.36
100	1.46	60	29.45
250	2.30		
500	3.26	(10) Unimetrics	
1ml	4.61mm	Series 9000	
2.5	7.28	10 μ l	0.46mm
5	10.30	25	0.73
10	14.57	50	1.03
25	23.03	100	1.46
50	32.57	250	2.30
		500	3.26
(5) Popper & Sons, Inc.		1000	4.61
Perfektum glass			
0.25ml	3.45mm		
0.5	3.45		
1	4.50		
2	8.92		
3	8.99		
5	11.70		
10	14.70		
20	19.58		
30	22.70		
50	29.00		

KDS Model 330 Series

Japanese Table of Syringe Diameters (Available in Japanese Models)

(1) "Air-Tite "All Plastic		(6) Nippro	
1cc	4.70mm	1ml long	4.75mm
2.5	9.70	1ml short	6.61
5.0	12.48	3	9.53
10	15.89	5	12.96
20	20.00	10	15.78
30	22.50	20	20.07
50	28.90	30	23.17
		50	29.13
(2) Becton Dickinson		(7) Hoshi	
Interim, WW design, Plastipak		1ml sm	4.80mm
1cc	4.70mm	1ml lg	6.70
3	8.59	2mm sm	6.70
5	11.99	2ml lg	9.20
10	14.48	3	10.30
20	19.05	5	12.20
30	21.59	10	15.00
60	26.60	20	19.00
		30	22.50
(3) Becton Dickson		50	25.50
Glass - all types		100	34.00
0.5cc	4.64mm	(8) Natsume	
1	4.64	0.25ml	2.60mm
2.5	8.66	0.50	3.20
5	11.86	1	4.30
10	14.34	2	6.30
20	19.13	3	7.30
30	22.70	5	9.50
60	28.60	(9) Top	
		1ml	4.70mm
(4) Hamilton		2	6.40
1000-Series Gastight		3	9.30
10µl	0.46mm	6	13.10
25	0.73	12	15.40
50	1.03	25	21.00
100	1.46	30	23.00
250	2.30	50	29.00
500	3.26	(10) Terumo	
1ml	4.61mm	1ml	4.73mm
2.5	7.28	3	9.00
5	10.30	5	13.04
10	14.57	10	15.79
25	23.03	20	20.18
50	32.57	30	23.36
		60	29.45
(5) JMC Air-Tite pls		(11) Terumo Japan	
1ml	4.66mm	1ml sm	4.73mm
2	6.90	1ml lg	6.50
2.5	9.10	3	8.95
5	12.62	5	13.00
10	14.34	10	15.80
20	19.68	20	20.15
30	22.44	30	23.10
50	28.80	50	29.10
100	36.68		

KDS Model 330 Series

Standard Minimum and Maximum Flow Rates

Syringe size	Diameter*	Minimum	Maximum
10 µl	0.46 mm	0.001 µl/hr	21.10 µl/min
25 µl	0.73 mm	0.003 µl/hr	53.15 µl/min
50 µl	1.03 mm	0.005 µl/hr	105.8 µl/min
100 µl	1.46 mm	0.009 µl/hr	212.6 µl/min
250 µl	2.3 mm	0.021 µl/hr	527.6 µl/min
500 µl	3.26 mm	0.042 µl/hr	1060 µl/min
1 ml	4.61 mm	0.083 µl/hr	2119 µl/min
2.5 ml	7.28 mm	0.207 µl/hr	5286 µl/min
3 ml	8.59 mm	0.288 µl/hr	7360 µl/min
5 ml	10.3 mm	0.414 µl/hr	634 ml/hr
10 ml	14.70 mm	0.842 µl/hr	1293 ml/hr
20 ml	19.05 mm	1.414 µl/hr	2171 ml/hr
30 ml	21.59 mm	1.817 µl/hr	2789 ml/hr
50 ml	28.9 mm	3.277 µl/hr	4998 ml/hr
60 ml	26.6 mm	2.757 µl/hr	4234 ml/hr
100 ml	34.9 mm	4.746 µl/hr	7289 ml/hr
140 ml	38.4 mm	5.746 µl/hr	8824 ml/hr

Syringes from different manufacturers can have slightly different limits.

Note: *This is a reference diameter used to calculate the flow rate. The specific diameter should be entered for your syringe type.

KDS Model 330 Series

Japanese Minimum and Maximum Flow Rates (Available in Japanese Models)

Syringe size	Diameter*	Minimum	Maximum
10 µl	0.46 mm	0.001 µl/hr	21.10 µl/min
25 µl	0.73 mm	0.003 µl/hr	53.15 µl/min
50 µl	1.03 mm	0.005 µl/hr	105.8 µl/min
100 µl	1.46 mm	0.009 µl/hr	212.6 µl/min
250 µl	2.3 mm	0.021 µl/hr	527.6 µl/min
500 µl	3.26 mm	0.042 µl/hr	1060 µl/min
1 ml	4.7 mm	0.087 µl/hr	2703 µl/min
2.5 ml	9.1 mm	0.323 µl/hr	6259 µl/min
3 ml	9.3 mm	0.337 µl/hr	8626 µl/min
5 ml	12.62 mm	0.7 µl/hr	953 ml/hr
10 ml	14.34 mm	0.802 µl/hr	1230 ml/hr
20 ml	19.68 mm	1.51 µl/hr	2317 ml/hr
30 ml	22.44 mm	1.963 µl/hr	3013 ml/hr
50 ml	28.8 mm	3.232 µl/hr	4963 ml/hr
60 ml	36.68 mm	5.243 µl/hr	8051 ml/hr

Syringes from different manufacturers can have slightly different limits.

Note: *This is a reference diameter used to calculate the flow rate. The specific diameter should be entered for your syringe type.

KDS Model 330 Series

Limited Warranty

KD Scientific Inc. warrants to the first consumer purchaser, for a period of one year from the date of purchase that this unit, when shipped in its original container, will be free from defective workmanship and materials and agree that it will, at its option, either repair or replace the defective unit.

This warranty does not extend to misuse, neglect or abuse, normal wear and tear, accident, modification or unauthorised repair.

KD Scientific will not be liable or in any way responsible for any incidental or consequential economic or property damage. Some States do not allow the exclusion of incidental or consequential damages, so the above exclusion may not apply to you.

There are no implied warranties of merchantability, or fitness for a particular use, or of any other nature. Some states do not allow this limitation on implied warranty, so the above limitation may not apply to you.

If a defect arises within the warranty period contact KD Scientific Inc., (see address below).

The customer is responsible for shipping charges and must first obtain a Return Material Authorization number (RMA) before the unit will be accepted. If a replacement unit is issued it is covered only for the remainder of the original warranty period dating from the purchase of the original device.

This warranty gives you specific legal rights. You may also have other rights which vary from state to state.

Note: This pump is not registered with the FDA and is not for clinical use on human beings.

Syringe pumps are manufactured by:

KD Scientific Inc.,

84 October Hill Road

Holliston, MA 01746

Phone: 508.429.6809 Fax: 508.893.0160

E-mail: info@kdsscientific.com

Web: www.kdsscientific.com

The logo for KD Scientific Inc. features the lowercase letters "kd" in a bold, blue, sans-serif font, followed by the word "Scientific" in a larger, blue, serif font. A horizontal line is positioned below the "kd" portion of the logo.

Copyright © 1997 by KD Scientific Inc. All rights reserved. Printed in USA

WEEE/RoHS Compliance Statement


EU Directives WEEE and RoHS

To Our Valued Customers:

We are committed to being a good corporate citizen. As part of that commitment, we strive to maintain an environmentally conscious manufacturing operation. The European Union (EU) has enacted two Directives, the first on product recycling (Waste Electrical and Electronic Equipment, WEEE) and the second limiting the use of certain substances (Restriction on the use of Hazardous Substances, RoHS). Over time, these Directives will be implemented in the national laws of each EU Member State.

Once the final national regulations have been put into place, recycling will be offered for our products which are within the scope of the WEEE Directive. Products falling under the scope of the WEEE Directive available for sale after August 13, 2005 will be identified with a "wheelie bin" symbol.

Two Categories of products covered by the WEEE Directive are currently exempt from the RoHS Directive - Category 8, medical devices (with the exception of implanted or infected products) and Category 9, monitoring and control instruments. Most of our products fall into either Category 8 or 9 and are currently exempt from the RoHS Directive. We will continue to monitor the application of the RoHS Directive to its products and will comply with any changes as they apply.



- **Do Not Dispose Product with Municipal Waste.**
- **Special Collection/Disposal Required.**