

To: «EDITOR2» «EDITOR3»

Publication: «PUBLICATION»



## NEW PRODUCT RELEASE

### **NEW Volumetric Infusion Syringe Pump features Simplistic Operation**

The **KD Scientific EZFlow 2040** is a **volumetric** infusion syringe pump designed to enhance quick efficient operation while maintaining simplicity. It meets the CE directive and is EC 93/42 EEC.

This pump is designed to deliver precise infusions at rates from 0.1 to 1200 ml/hr. It features a bright 2.7 inch OLED with a separate display which shows volume delivered, volume rate, time and date.

There are 5 preset languages and many other advanced features that make the **KDS EZFlow 2040** one of the most outstanding infusion pumps in its class.

The **KDS EZFlow 2040** is portable, user friendly and offers a new alternative for micro reliable infusions. The ergonomic touch screen panel provides ease of selection and programming.

The pump has a self-diagnostic program which initiates each time the pump is powered on. Once completed successfully the pump is ready for use.

The **KDS EZFlow 2040** has an AC power source with battery back up of approximately 3 hours continuous operation at 25 ml/hr. It has a built in rechargeable Ni-MH battery and a RS232 serial port for data communication. The pump can be used on a horizontal platform or optional pole mount adapter for maximum flexibility.

**KD Scientific** designs, manufactures and sells a range of quality fluidics equipment used by research laboratory markets worldwide.



KD Scientific syringe pumps are an economical solution to delivering precise and smooth flow in research, pilot plants and production applications. They are recognized worldwide for quality, accuracy and reliability. A broad line of syringe pumps are offered: from a simple one syringe infuse only, to a programmable multi-syringe infuse/withdrawal pump. **KD Scientific** operates as a separate business within Harvard Bioscience at their Holliston, Massachusetts facility.

**For more information contact:**

**KD Scientific**

84 October Hill Road  
Holliston, Massachusetts 01746

Phone: 508-429-6809

Fax: 508-893-0160

E-mail: [info@kdscientific.com](mailto:info@kdscientific.com)

Web: [www.kdscientific.com](http://www.kdscientific.com)

Product and company names listed are trademarks or trade names of their respective companies.

**NOTE:** This Press Release has been sent electronically. If it has been sent to the wrong individual, or if you no longer wish to receive KD Scientific Press Releases, please notify us at: [PressRelease@optonline.net](mailto:PressRelease@optonline.net)