



NEW PRODUCT RELEASE

Amazing *NEW* Spill Sensor Detects Leaks

KD Scientific has released the **NEW OS 250**, a system which detects spills and leaks before they cause a problem.

As little as 3 drops of liquid will cause the OS 250 to react.

The system consists of a moisture sensing mat, a material specially developed for detecting liquid spills. A simple connection between the mat and the OS 250 Controller will sound an audible alarm, flash an LED and will turn off the power of any device plugged into the single outlet, solid state power controller. The switched power outlet can control up to 8 amps.



The **OS 250 Spill Sensor** is supplied with the controller and four reusable 30 x 30 cm mats that can be cut to any size with a sharp scissors or knife. It also includes the connector cable between the mat and the controller. The overall dimensions of the controller are 4.6 x 8.2 x 18 cm.

Applications for the OS 250 are in pharmaceutical, chemical, petrochemical, biotechnology, semiconductor, plastics, industrial, government, and scientific research and development markets.

KD Scientific Inc. designs, manufactures and sells a wide range of quality fluidics equipment used by research laboratory markets worldwide.

For more information contact:

KD Scientific Inc.

84 October Hill Road
Holliston , Massachusetts 01746

Phone: 508-429-6809

Fax: 508-893-0160

E-mail: info@kdscientific.com

Web: www.kdscientific.com