

To: «EDITOR2» «EDITOR3»
Publication: «PUBLICATION»



NEW PRODUCT RELEASE

Protect your Microscope from Spills!

KD Scientific has introduced the *NEW OS-250*, a system which detects spills and leaks before they cause a problem in an expensive microscope or other piece of equipment.

As little as 3 Drops of Liquid will cause the **OS-250** to react. The system consists of a moisture sensing mat and control unit. The mat is made from a material specially developed for detecting liquid spills. It is connected to the control unit by a simple connection cable.

When liquid is detected on the sensing mat the **OS-250 Controller** will sound an audible alarm, flash an LED and will turn off the power of any device plugged into the single outlet, solid state power controller. The switched power outlet can control up to 8 amps.

The **OS-250 Spill Sensor** is supplied with the controller and four reusable 30 x 30 cm mats that can be cut to any size with a sharp scissors or knife. It also includes the connector cable between the mat and the controller.

Applications for the **OS-250 Spill Detector** exist in any areas where spills will cause problems or damage to equipment.

KD Scientific designs, manufactures and sells a wide range of quality fluidics equipment used by research laboratory markets worldwide.

For more information contact:

KD Scientific

84 October Hill Road
Holliston, Massachusetts 01746

Phone: 508-429-6809

Fax: 508-893-0160

E-mail: info@kdsscientific.com

Web: www.kdsscientific.com



Product and company names listed are trademarks or trade names of their respective companies.

NOTE: This Press Release has been sent electronically. If it has been sent to the wrong individual, or if you no longer wish to receive KD Scientific Press Releases, please notify us at: PressRelease@optonline.net